

Nurse Console Soft Key Command Guide

The following guide identifies all soft key functions available on the Nurse Console. Not all functions may be enabled for your system. Soft Key label descriptions and colors may also vary from the defaults shown below. Soft Keys may be toggled On/Off, or can be used to initiate a specific function such as Set Time.

All Page	An audio page is sent to all rooms in the selected coverage area for the console.
All Page Site	An audio page is sent to all rooms in the entire facility.
All Page Urgent	An audio page is sent to all rooms in the selected coverage area and any existing audio connections are preempted.
Bed Attn	Upgrades a bed station pillow speaker button's to the next highest priority level.
Bed Patient	Downgrades a bed station pillow speaker button's back to the previous priority level.
Bed Special	Upgrades a bed station pillow speaker button's to the Special status.
BedUp Review	Reviews the bed upgrade status for all beds in the coverage area.
EVS Staff	Identifies all rooms where EVS staff members are registered.

Follow	Enables call tones to sound in all patient rooms with registered staff members.
Home	Returns to the main soft-key screen.
LPN Staff	Identifies rooms with staff members registered as LPN staff level.
Monitor	Pressing this key allows entry of room numbers for audio monitoring. Format is Area+Room Number:Bed Number. Use a comma to separate multiple rooms. Rooms set to Privacy cannot be monitored.
Mute Tones	The active call tone is muted until another call becomes active. Normal tone volume is restored after 2 minutes if no other calls become active.
Need EVS	Sets a Service reminder for EVS staff level for the selected room number.
Need LPN	Sets a Service reminder for the LPN staff level for the selected room number.
Need PCT	Sets a Service reminder for the PCT staff level for the selected room number.

**Need
RN**

Sets a Service reminder for the RN staff level for the selected room number.

Night

Toggles On/Off the Night tone level for tone annunciation only at the selected console. Day Tone level is restored when toggled Off.

**PCT
Staff**

Identifies rooms with staff members registered as PCT staff level.

Privacy

Sets a room into Privacy mode to prevent monitoring room audio sounds. Format is Area+Room Number:Bed Number. Use a comma to separate multiple rooms.

**RN
Staff**

Identifies rooms with staff members registered as RN staff level.

SIP

Enables or disables SIP phone extensions. Calls are not directed to SIP phones that are disabled.

**Site
Night +**

Enables the night tone volume for all consoles, annunciator panels, and duty stations within the entire facility.

**Site
Night -**

Disables the night tone volume for all consoles, annunciator panels, and duty stations within the entire facility.

**Staff
Page**

An audio page is sent to all rooms in the selected coverage area with registered staff members.

**Staff
Page
Urgent**

An audio page is sent to all rooms in the selected coverage area with registered staff members and any existing audio connections are preempted.

Time

Sets system time for the nurse call system including time displayed on consoles and Reports.

Urgent

Upgrades the selected call to Urgent status, escalating the call to a higher escalation chain.

Vol +

Increases the volume of the room station which is in communication with the console.

Vol -

Decreases the volume of the room station which is in communication with the console.

Refer to technical manual KI-2907 *Responder® 5000 Component Installation Guide* for additional information including all warnings and cautions regarding console use.