



Assembly Test Form for Monitor Calibration - Service

1 Test Execution Summary

Test Date:	6/07/2023
Tester (print & sign):	JORGE F. + DENIS C.
Monitor S/N:	6127 2220 11

2 Test Equipment Used

Test Equipment	Model	Part Number	Serial Number
EIZO Monitor Scripts	N/A	6456	N/A
USB 3.0 Type A Male to B SuperSpeed Cable	N/A	N/A	N/A
EIZO UX2 Sensor	N/A	7042	N/A SAV/B67 2242716
Colorimeter	X-Rite	I1 displaypro	N/A
USB Keyboard	N/A	N/A	N/A

Notes:

- All test equipment shall have up-to-date calibration in accordance with ISO 13485:2016 requirements.
- Observe all ESD rules as per SOP Electrostatic Discharge Controls. Check wrist strap before starting.

3 Test Results

Test number corresponds to Assembly Test Procedure number.

3.1 ExactVu Monitor Calibration

Test #	Pass Criteria	Pass/ Fail	Comments
5.1	Check Monitor Date of Manufacture		Mfd: <u>2021-05-31</u> If 2018-12 or later, run sect. 5.1.1 If 2018-11 or earlier, run sect. 5.1.2
5.1.1.1 ONLY for monitors DOM 2018-12 or later	Choose one of the following: <input type="checkbox"/> RadiCS software already installed <input checked="" type="checkbox"/> RadiCS software installation required; Installation performed.	P	
	UX2 Sensor plugged in and detected	P	
	Checked System Hardware Version		Please circle one: 1) 2.6 or 2.9 2) 2.5
	USB 3.0 Type A Male to B SuperSpeed Cable plugged in to the monitor and system. Skip if the system has 2.6 or 2.9 HW		
	Acceptance Test Passed	P	Note: If test fails, section 5.1.1.2 must be run
	This section ONLY needs to be run if the above acceptance FAILS or a new monitor was installed. And ONLY for monitors DOM 2018-12 or later		
5.1.1.2 ONLY for monitors DOM 2018-12 or later	USB 3.0 Type A Male to B SuperSpeed Cable plugged in to the monitor and system. Skip if the system has 2.9 HW	P	
	"AllLockOff.bat" script executed	P	
	Monitor Menu -> Color and set: • Mode to "User1" • Brightness to "70" • Contrast to "50" • Temperature to "7500k" • Gamma to "2.2"	P	
	Go to Control Panel -> Appearance and Personalization -> Display -> Change Display Settings -> Advanced Settings -> Color Management -> Color Management ... All Profiles are removed "Use my settings for this device" is unchecked	P	
	Calibration Target set as: • Lmax: "160.00" • Lmin: check "Set Lmin as low as possible" • Color: "7500" • Display Function: "DICOM Part 14 GSDF"	P	
	QC Guideline set to "Basic QC Secondary"	P	

	UX2 Sensor plugged in and detected	P	
	Close the RadiCS tool	P	
	Monitor Menu -> Color, record settings		Mode: <i>V21</i> Brightness: <i>94</i> Contrast: <i>50</i> Temperature: <i>7500K</i> Gamma: <i>22</i>
	Reboot the system	P	
	Exit to the Shell and kill the ExactVu App	P	
	Monitor Menu -> Color, set the below to the values recorded previously: <ul style="list-style-type: none"> • Mode: • Brightness: • Contrast: • Temperature: • Gamma: 	P P P P P	
	"AllLockOn.bat" script executed	P	
	All monitor buttons are locked including the power button	P	
	Acceptance Test Passed	P	
5.1.2.1 ONLY for monitors DOM 2018-11 or earlier	Choose one of the following: <input type="checkbox"/> i1Profiler already installed <input type="checkbox"/> i1Profiler installation is required. Installation performed.	NA	
5.1.2.2 ONLY for monitors DOM 2018-11 or earlier	I1 Display Pro Colorimeter plugged in. i1Profiler Software started and meter detected	NA	
	ExactVuDICOM.dpxf loaded. Settings verified.	NA	
	Colorimeter aligned and placed into the center of the monitor	NA	
	Measurements completed. ExactVuDICOM.icm profile created	NA	
	Results within tolerances: <ul style="list-style-type: none"> • White Point: 7500K +/- 200 • Luminance: 160 cd/m² +/- 5 	NA	White Point result: _____ Black Luminance result: _____ Luminance result: _____
	Quality Measurement passed	NA	

D. CHARLET

