



report

NO. H3192

INFORMATION FOR DISTRIBUTORS

August 11, 2003

TO: ALL HEALTH CARE DISTRIBUTORS – U.S., CANADA & EXPORT

**SUBJECT: CALL CORD MODELS CC200, CY210, AND CCDIN UL
APPROVED FOR USE IN OXYGEN RICH ENVIRONMENT**

The following call cords for the Responder® Models CC200, CY210, and CCDIN have been tested and approved for use in oxygen rich environments by Underwriter Laboratories in accordance with the Hospital Signaling and Nurse Call Equipment standard (UL 1069). The approved call cord models now carry the label - "Note: May be used by patients undergoing oxygen therapy – hang on hook (in holder) when not in use."



An important concern in hospitals today is patient safety. Rauland-Borg has the unique advantage of offering a call cord that can be used by a patient receiving oxygen treatment. From a risk management perspective, the hospital has the reassurance that their patients receiving oxygen treatment can safely use the pillow speakers to control the television, lights, and call and speak to the nursing staff. No other nurse call manufacturers can offer the same reassurance. Rauland-Borg is the ONLY manufacturer of Health Care communications products to obtain this UL approval.

RAULAND-BORG CORPORATION

3450 West Oakton Street, Skokie, IL 60076-2951 • Tel: (847) 679-0900 • FAX (847) 679-4106

In Canada: RAULAND-BORG (CANADA) Inc. • 4025 Sladeview Crescent, Units 4-6, Mississauga, ON, Canada L5L 5Y1 • (905) 607-2335 • FAX: (905) 607-3554