

Overview

Computer name: TRINITY-363
Software Version: 4.3.3 (R1)
IK: 9CF3A785E4059594
Security Code: E4059594
Partition C: Free: 174.457Go Total: 199.512Go
Partition D: Free: 1618.5Go Total: 1663.01Go
Host name: TRINITY-363
IP: 127.0.0.1

\$SYSTEM_CONFIG_TITLE\$

\$SYSTEM_CONFIG_CONTENT\$

Application Config

APPLICATION_MIRROR: 0
AllowBiopsy: 1
AllowFastRegistration: 1
AllowFly: 0
AllowMeasureCalculation: 1
AllowPatientExport: 0
AllowTreatment: 0
AutoValidationTimeout: 0
IK: 9CF3A785E4059594
KeepMrExamInReceptionDirectoryDelay.Timeout: 3
ListLocalFolderPathsForExport: 0
ListNetworkPathsForExport: 0
MaxLinesInCommentEditor: 20
PatientManagerRootDir: D:\KOELIS
SaveUSNativeVolume: 0
TransPerinealSupported: 0
TransRectalSupported: 1
VirtualKeyboardActivation: 1

Users Config

DEMO

2ndLookOptionSupported: 1
ActivatedGuides: 5 KRNG.EL1 KRNG.EL2 KRNG.EL4 KSU.EL1 KSU.EL4
AnyModalitiesFusionWidgetSupported: 1
ApplicationLegacyDiagnosis: 0

ApplicatonHardShutdownPreviousTime: 0
DataSegmentationDialog.LinkMR: 0
Demo: DemoDataSet
ExportOption_SC_CARTOGRAPHY: 1
ExportOption_SC_PDFREPORT: 1
ExportOption_SC_SCREENSHOTS: 1
ExportOption_SC_SCREENSHOTS_UNCOMPRESSED: 0
ExportOption_SURFACESEGMENTATION: 0
IdClient: DEMO
ImportOption_CT: 0
ImportOption_MR: 1
ImportOption_RTSTRUCT: 1
ImportOption_SECONDARYCAPTURE: 0
ImportOption_SEGM: 1
ImportOption_TEP: 0
ImportOption_USENHANCED: 1
Lang: es_es
LastFirstLookVisitedExam: 2602\20220328-002656-TESTAM\202203280026\
LogOffTimeout: 60
LogOffTimeoutMax: 240
LogOffTimeoutMin: 10
NavigationOptionDisplayBox.InitialValues.Show2ndLookData: 1
NavigationOptionDisplayBox.InitialValues.ShowCurrentExam: 548
NavigationOptionDisplayBox.InitialValues.ShowGenerals: 3342336
NavigationOptionDisplayBox.InitialValues.ShowPlanningExam: 802
NavigationOptionDisplayBox.InitialValues.ShowPreviousExam: 545
PanoramaAcquisitionMode: PANORAMA_TRANSRECTAL_3
PeropImageViewerVisuLeft.ChangeAutomatically: 0
PlanningGuide2D: 0
PlanningOtherGuide: 0
TrusMrFusionSupported: 1
VirtualKeyboard.DefaultLanguage: es_es
VirtualKeyboard.SupportedLanguages: 1 es_es

PROMAP

2ndLookOptionSupported: 1
ActivatedGuides: 5 KRNG.EL1 KRNG.EL2 KRNG.EL4 KSU.EL1 KSU.EL4
AnyModalitiesFusionWidgetSupported: 1
ApplicationLegacyDiagnosis: 0
ApplicatonHardShutdownPreviousTime: 0
DataSegmentationDialog.LinkMR: 0
Demo:
ExportOption_SC_CARTOGRAPHY: 1
ExportOption_SC_PDFREPORT: 1
ExportOption_SC_SCREENSHOTS: 1
ExportOption_SC_SCREENSHOTS_UNCOMPRESSED: 0
ExportOption_SURFACESEGMENTATION: 0
IdClient: PROMAP
ImportOption_CT: 0
ImportOption_MR: 1
ImportOption_RTSTRUCT: 1
ImportOption_SECONDARYCAPTURE: 0
ImportOption_SEGM: 1

ImportOption_TEP: 0
ImportOption_USENHANCED: 1
Lang: es_es
LastFirstLookVisitedExam: 2602\20230313-105510-COÑA H\202306081508\
LogOffTimeout: 60
LogOffTimeoutMax: 240
LogOffTimeoutMin: 10
NavigationOptionDisplayBox.InitialValues.Show2ndLookData: 1
NavigationOptionDisplayBox.InitialValues.ShowCurrentExam: 548
NavigationOptionDisplayBox.InitialValues.ShowGenerals: 3342336
NavigationOptionDisplayBox.InitialValues.ShowPlanningExam: 802
NavigationOptionDisplayBox.InitialValues.ShowPreviousExam: 545
PanoramaAcquisitionMode: PANORAMA_TRANSRECTAL_3
PeropImageViewerVisuLeft.ChangeAutomatically: 0
PlanningGuide2D: 0
PlanningOtherGuide: 0
TrusMrFusionSupported: 1
VirtualKeyboard.DefaultLanguage: es_es
VirtualKeyboard.SupportedLanguages: 1 es_es

Reports Config

DEMO

DEPARTMENT: 1 Urology
HOSPITAL: 5 CENCOMEX S.A Galvarino 7640, Parque Aconcagua- Quilicura Santiago Chile
LOGO: cencomex.jpg

PROMAP

DEPARTMENT: 1 Urology
HOSPITAL: 5 CENCOMEX S.A Galvarino 7640, Parque Aconcagua- Quilicura Santiago Chile
LOGO: cencomex.jpg

Service class providers

MR

AETitle: KOELIS
Cli: -su S --fork
ConfigurationFile: Storescp.cfg
ConfigurationProfile: MR_ONLY
Name: MR
Port: 11112
ReceiptDirectory: D:\temp\KOELIS\0000-KoelisDicom\MR



Service class users