

Overview

Computer name: TRINITY-168
Software Version: 4.5.0 (R3)
Revision Id: 15fc5a479d1fd39a349924ff89de945ea3e202ca
Security Code: 1ADF-10FB
Partition C:\: Free: 159.514Go Total: 195.312Go
Partition D:\: Free: 1628.37Go Total: 1667.48Go
Partition F:\: Free: 1.849Go Total: 1.86218Go
Application installation path : "C:\Windows\system32"
Version of configuration files: 4.5.0 (R1)
Host name: TRINITY-168
IP: 127.0.0.1

Application Config

APPLICATION_MIRROR: 0
AutoValidationTimeout: 7
CounterTempPeriods: 3 1 0 0
IsTrinity: 1
KeepMrExamInReceptionDirectoryDelay.Timeout: 2
ListLocalFolderPathsForExport: 0
ListNetworkPathsForExport: 0
MaxLinesInCommentEditor: 20
PatientManagerRootDir: D:\KOELIS
SaveUSNativeVolume: 0
UltrasoundScannerSupported: 2 Medison Samsung
VirtualKeyboardActivation: 1
BIOPSY: 1
PROMAP-IC: 0
PROMAP-FT: 0
FLY: 0
FAST: 1
PROMAP-VM: 0
PROMAP-TP: 0
PROMAP-TR: 1
PROMAP-EX: 0

Users Config

DEMO

ApplicationLegacyDiagnosis: 0
ApplicatonHardShutdownPreviousTime: 0
DataSegmentationDialog.LinkMR: 0
Demo: DemoDataSet
ExportOption_SURFACESEGMENTATION: 0

IdClient: DEMO
ImportOption_SECONDARYCAPTURE: 0
Lang: es_es
LastFirstLookVisitedExam:
LogOffTimeout: 60
LogOffTimeoutMax: 240
LogOffTimeoutMin: 10
NavigationOptionDisplayBox.InitialValues.Show2ndLookData: 1
NavigationOptionDisplayBox.InitialValues.ShowBiopsies: 1
NavigationOptionDisplayBox.InitialValues.ShowBiopsiesNearAxialPlane: 0
NavigationOptionDisplayBox.InitialValues.ShowCancerousBiopsiesOnly: 0
NavigationOptionDisplayBox.InitialValues.ShowCurrentExam: 548
NavigationOptionDisplayBox.InitialValues.ShowGenerals: 3342336
NavigationOptionDisplayBox.InitialValues.ShowGrid: 0
NavigationOptionDisplayBox.InitialValues.ShowImages: 1
NavigationOptionDisplayBox.InitialValues.ShowPlanningExam: 802
NavigationOptionDisplayBox.InitialValues.ShowPreviousExam: 545
NavigationOptionDisplayBox.InitialValues.ShowProstate: 1
NavigationOptionDisplayBox.InitialValues.ShowVirtualBiopsyInRake: 1
PanoramaAcquisitionMode: PANORAMA_TRANSRECTAL_3
PeropImageViewerVisuLeft.ChangeAutomatically: 0
UltrasoundImage.ProbeDirectionFlip: 0
UltrasoundScanner: Medison
VirtualKeyboard.DefaultLanguage: es_es
VirtualKeyboard.SupportedLanguages: 1 es_es

PROMAP-2L: 1
ELASTIC FUSION: 1
REGISTRATION WIDGET: 1
PROMAP-GR: 0
PROMAP-PL: 0
MR: 1
US UNHANCED: 1
CT: 1 validity period: [6/23/2023,7/14/2023]
TEP: 1 validity period: [6/23/2023,7/14/2023]
RTSTRUCT (KOELIS): 1
SEGM (KOELIS): 1
KRNK.EL1: 1
KRNK.EL2: 1
KRNK.EL4: 0
KSU.EL1: 0
KSU.EL4: 0
KRG.EL5: 0
3D MAPPING (DICOM): 1
REPORT (DICOM): 1
SCREENSHOT (DICOM): 1
SCREENSHOT (DICOM) Uncompressed: 0

PROMAP

ApplicationLegacyDiagnosis: 0
ApplicatonHardShutdownPreviousTime: 0
DataSegmentationDialog.LinkMR: 0
Demo:
ExportOption_SURFACESEGMENTATION: 0

IdClient: PROMAP
ImportOption_SECONDARYCAPTURE: 0
Lang: es_es
LastFirstLookVisitedExam: 2602\20230613-172128-PEREZ\202306131721\
LogOffTimeout: 60
LogOffTimeoutMax: 240
LogOffTimeoutMin: 10
NavigationOptionDisplayBox.InitialValues.Show2ndLookData: 1
NavigationOptionDisplayBox.InitialValues.ShowBiopsies: 1
NavigationOptionDisplayBox.InitialValues.ShowBiopsiesNearAxialPlane: 0
NavigationOptionDisplayBox.InitialValues.ShowCancerousBiopsiesOnly: 0
NavigationOptionDisplayBox.InitialValues.ShowCurrentExam: 548
NavigationOptionDisplayBox.InitialValues.ShowGenerals: 3342336
NavigationOptionDisplayBox.InitialValues.ShowGrid: 0
NavigationOptionDisplayBox.InitialValues.ShowImages: 1
NavigationOptionDisplayBox.InitialValues.ShowPlanningExam: 802
NavigationOptionDisplayBox.InitialValues.ShowPreviousExam: 545
NavigationOptionDisplayBox.InitialValues.ShowProstate: 1
NavigationOptionDisplayBox.InitialValues.ShowVirtualBiopsyInRake: 1
PanoramaAcquisitionMode: PANORAMA_TRANSRECTAL_3
PeropImageViewerVisuLeft.ChangeAutomatically: 0
UltrasoundImage.ProbeDirectionFlip: 0
UltrasoundScanner: Medison
VirtualKeyboard.DefaultLanguage: es_es
VirtualKeyboard.SupportedLanguages: 1 es_es

PROMAP-2L: 1
ELASTIC FUSION: 1
REGISTRATION WIDGET: 1
PROMAP-GR: 0
PROMAP-PL: 0
MR: 1
US UNHANCED: 1
CT: 1 validity period: [6/23/2023,7/14/2023]
TEP: 1 validity period: [6/23/2023,7/14/2023]
RTSTRUCT (KOELIS): 1
SEGM (KOELIS): 1
KRNGL1: 1
KRNGL2: 1
KRNGL4: 0
KSUL1: 0
KSUL4: 0
KRGAL5: 0
3D MAPPING (DICOM): 1
REPORT (DICOM): 1
SCREENSHOT (DICOM): 1
SCREENSHOT (DICOM) Uncompressed: 0

Reports Config

DEMO

DEPARTMENT: 1 Urologia

HOSPITAL: 4 Munich Pharma Medical LTDA ZI Salto 4875 Huechuraba 7810808 Santiago
CHILE
LOGO: LOGO-SOCO FAR.jpg

PROMAP

DEPARTMENT: 1 UROLOGIA
HOSPITAL: 4 UROFUSION La Concepcion 191, of. 201 Telefono: 993326067 Email:
contacto@urofusion.cl
LOGO: Logo Urofusión.PNG

Service class providers

MR

AETitle: KOELIS
Cli: -su S --fork
ConfigurationFile: Storescp.cfg
ConfigurationProfile: MR_ONLY
Name: MR
Port: 11112
ReceiptDirectory: D:\temp\KOELIS\0000-KoelisDicom\MR

Service class users

AETitle: KOELIS
Cli: --no-halt +sd +r -q
ConfigurationFile: Storescu.cfg
ConfigurationProfile: MR_ONLY
IpAddress: 192.168.109.200
IsActive: 1
Name:
Port: 11112
Role.Query_Retrieve: 0
Role.Storage: 1
Role.Worklist: 0