

Please fill in the first start up form and send it back to Quanta System Service Dept.via e-mail at [service@quantasystem.com](mailto:service@quantasystem.com) (scanned version)

### FIRST START-UP FORM

Model of the laser system:

☐ **Cyber Ho 100**      ☐ **Litho 100**      ☐ **Cyber Ho 150**      ☐ **Litho 150**

Serial number: \_\_\_\_\_

Nominal Line Voltage \_\_\_\_\_ V, \_\_\_\_\_ Hz

<b>Checks</b>	<b>Put x if performed</b>	<b>Comments</b>
Visual check of the laser system		
Integrity and presence of the Accessories		
Integrity and presence of the Power cord		
Functioning of the Touchscreen		
Fill the Cooling system		
Presence of the Aiming beam out of the delivery system		
Absence of warnings or alarms		

General comments after the installation: \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_

Distributor Company Name: \_\_\_\_\_

Country: \_\_\_\_\_

Check performed by: \_\_\_\_\_

Date / Signature: \_\_\_\_\_