

Please fill in the following first start up form and send it back to Quanta System Service Dept.:
via e-mail at service@quantasystem.com (scanned version)

FIRST START-UP FORM

Model of the laser system: **Cyber Ho 60/ Litho 60**

Serial number: _____

Nominal Line Voltage _____ V, _____ Hz

Checks	<i>Put x if performed</i>	Comments
Visual check of the laser system		
Integrity and presence of the Accessories		
Integrity and presence of the Power cord		
Functioning of the Touchscreen		
Fill the Cooling system		
Presence of the Aiming beam out of the delivery system		
Absence of warnings or alarms		

General comments after the installation: _____

Distributor Company Name: _____

Country: _____

Check performed by: _____

Date / Signature: _____