

3.8 First start-up form

Please fill in the following first start-up form and send it back via e-mail to service@quantasystem.com (scanned version) or fax it to +39 0331 367815.

FIRST START-UP FORM

Model of the laser system: **Litho DK30**

Serial number: _____

Line Voltage: 230V 110V (please circle the correct voltage)

Customer's name and address: _____

Contact phone number/email address: _____

Checks	Put x if performed	Comments
Visual check of the laser system		
Accessories		
Main ground wires		
Electrical connections		
Touchscreen		
Cooling system		
Aiming beam		
The laser beam alignment		
Energy check		

General comments after the installation: _____

Authorized company for technical service: _____

Country: _____

Service technician/engineer: _____

Signature: _____

Please send the scanned first start-up form to service@quantasystem.com or fax it to +39 0331 367815.