

Please fill in the following Installation Report and send it back to Quanta System Service Dept.:

via e-mail at service@quantasystem.com (scanned version),
or fax it to +39 0331 367815.

FIRST START-UP FORM

Model of the laser system: **Fiber Dust**

Serial number: _____

Nominal Line Voltage _____ V, _____ Hz

<i>Checks</i>	<i>Put x if performed</i>	<i>Comments</i>
Visual check of the laser system		
Integrity and presence of the Accessories		
Integrity and presence of the Power cord		
Functioning of the Touchscreen		
Presence of the Aiming beam out of the delivery system		
Absence of warnings or alarms		

General comments after the installation: _____

Distributor Company Name: _____

Country: _____

Check performed by: _____

Date / Signature: _____