


# Qualitas Series Parts List

This parts list contains the parts that are required for repairing and maintaining Qualitas Series.

- Please note that the listed parts are the smallest unit that can be provided.
- Please read the instruction manual to double check the parts name and parts number.
- Due to improvement of the product, parts of this parts list may differ slightly from the product purchased.
- If you have any question, please contact your dealer.



PARAMOUNT BED CO., LTD.

## Contents

	Page
1. Parts List for Qualitas $\alpha$ (3 MOTOR TYPE) . . . . .	3
2. Parts List for Qualitas $\beta$ (4 MOTOR TYPE) . . . . .	4

## How to use parts list

- The color or dimension differences do not reflect on the diagrams.  
Please check the part table for correct unit color and size.
- Each part is numbered. The number in picture and the number in chart correspond to each other.
- Be sure to inform us the parts number for parts inquiry.

## When ordering repair parts

- Each parts will be ordered by unit of one.
- When ordering, please inform us which part has a problem, the parts number and quantity required.

## Repair during the warranty period

While the most repairs are carried out free of charge as specified in the sales contract, some types of repairs are carried out at the user's expense even during the warranty period. See the sales contract for more information.

## After the warranty

When problems can be solved by repairs, they can be carried out at the user's expense.

## Minimum part availability period

Paramount Bed will stock bed replacement parts for repair (needed to maintain product performance) for at least 6 years after production has been discontinued.

## Cautions

- For user's convenience, Paramount Bed try to provide the parts in details as much as possible. However, due to safety and quality aspect and efficiency for delivery, some parts may be delivered in assembled form.
- Paramount Bed do not provide packing materials. (cardboard, tape, band and etc) Please prepare them at customer's cost.
- Paramount Bed try to deliver the replacement parts quickly after your order, but sometimes it might take a time for delivery due to distribution conditions. Please check the delivery time when ordering.
- When you replace the parts, Paramount Bed recommend to replace fixing parts together with the replacing parts to prevent looseness.  
(some parts are not suitable for reuse)
- When replacing the parts, please make sure they are fixed by the screws, bolts and nuts.





### Caution

- When fixing, please do not apply too much pressure on the screws and bolts. Too much tightening may cause deformation or damage.
- Be careful handling snap-pins. When you attach or detach, there is a possibility that snap-pin may pop out in unexpected direction and cause an injury.

## Parts List for Hospital Electrical Bed Qualitas $\alpha$ (3M)

No.	Product code	Description
01	3A10000500A	Back section base
02	3A10001300A	Hip section base
03	3A10001800A	Knee section base
04	3A10001900A	Leg section base
05	3F01237900A	Mattress stopper (Mattress Length 1,910mm)
	3F07300300A	Mattress stopper (Mattress Length 2,000mm)
06	3A50004300A	H board AS (Blue)
	3A50008300A	H board AS (Medium Oak)
	3A50008500A	H board AS (Green)
07	3A50004400A	F board AS (Blue)
	3A50008400A	F board AS (Medium Oak)
	3A50008200A	F board AS (Green)
08	3F06331200A	Corner bumper
09	3A50006901A	Bumper attachment (Right)
10	3A50006902A	Bumper attachment (Left)
11	3C01115801A	Side rail (Right)
12	3C01115802A	Side rail (Left)
13	T125LNS	Caster
14	T125LMS	Caster (Steering)
	T125LMSX	Anti-static caster (Steering) *Option
15	3F06824100A	Caster lock pedal (Right)
16	3F06824000A	Caster lock pedal (Left)
17	TA10000500A	Actuator (Back raise)
18	TA10000300A	Actuator (Knee raise)
19	TA10000600A	Actuator (Height up)
20	TC10000100A	Control Box with cable power supply
21	TH20000200A	Hand controller



## Parts List for Hospital Electrical Bed Qualitas $\beta$ (4 Motor)

No.	Product code	Description
01	3A10000500A	Back section base
02	3A10001300A	Hip section base
03	3A10001800A	Knee section base
04	3A10001900A	Leg section base
05	3F01237900A	Mattress stopper (Mattress Length 1,910mm)
	3F07300300A	Mattress stopper (Mattress Length 2,000mm)
06	3A50004300A	H board AS (Blue)
	3A50008300A	H board AS (Medium Oak)
	3A50008500A	H board AS (Green)
07	3A50004400A	F board AS (Blue)
	3A50008400A	F board AS (Medium Oak)
	3A50008200A	F board AS (Green)
08	3F06331200A	Corner bumper
09	3A50006901A	Bumper attachment (Right)
10	3A50006902A	Bumper attachment (Left)
11	3C01115801A	Side rail (Right)
12	3C01115802A	Side rail (Left)
13	T125LNS	Caster
14	T125LMS	Caster (Steering)
	T125LMSX	Anti-static caster (Steering)
15	3F06824100A	Caster lock pedal (Right)
16	3F06824000A	Caster lock pedal (Left)
17	TA10000700A	Actuator (Back raise)
18	TA10000800A	Actuator (Knee raise)
19	TA10000400A	Actuator (Height up)
20	TC10000200A	Control Box with cable power supply
21	TH20000300A	Hand controller
22	TNP10001	Nurse control panel
23	3A100189	Battery Unit

