

Performance Test Litho

Model of the system: LITHO


Serial Number: LHT 0327 - 0219

Optical Test (with user autoadjust on)

Fiber	Hz / J	Power on the system display	Power read on the power meter	Lamp (us)
Lithotripsy mode				
365	10Hz / 3.0J	30 W	28.8 W	631
365	10Hz / 1.2J	12 W	12.5 W	418
365	12Hz / 1.0J	12 W	12.1 W	382
200	10Hz / 1.2J	12 W	12.5 W	418
200	12Hz / 1.0J	12 W	12.4 W	382

If the results are in max 5% different you can do the next step

Long run test (10 minutes lasering in continuous Lithotripsy mode)

200	12Hz / 1.0J	Temperature of the fiber jack	cold ☒	warm ☐	hot ☐
Tested By: Jonge I	Surname: FERNANDEZ	Signature: 	Date: 19/5/2023		